


| Please check one: | | | | | | | | |
|-------------------|-------|--|--|--|--|--|--|--|
| Quotation | Order | | | | | | | |

Custom Quotation/Order for Solvent Desorption Sorbent Tubes

Please complete each numbered section on both pages of this specification form and send to your SKC representative:

| | SKC Inc. Tel: 800-752-8472 skcorder@skcinc.com | SKC South Tel: 800-752-7684 skcorder@skcinc.com | SKC-West Tel: 800-752-9378 custserv@skcwest.com | SKC Limited Tel: (01258) 480188 info@skcltd.com | SKC South Africa Tel: 11 913 2666 sales@skcsouthafrica.co.za | SKC Asia Tel: (65) 6271 0291 sales@skc-asia.com |
|-------|--|---|---|---|--|---|
| | pe of Order: (Please | | a tula a | | | |
| _ | New order for so Repeat order (P | - | | |) | |
| 2. Qu | antity ordered | (Note: | All custom orders | are subject to a : | ± 10% variation in qua | antity) |
| 3. An | alyte | 4. / | Analytical method | l number and ag | ency (if available) | |
| 5. Tu | be specifications (| Mark any specia | l instructions on o | drawing below, i | ncluding the blank sp | pace in the tube) |
| A. | Glass: | | | | | |
| | O.D.: Select One | | Other (| please specify): | | |
| | Length: Select Or | ne | Other (| please specify): | | |
| | Sealed tip: Selec | et One | | | | |
| | Additional Notes | : : | | | | |
| | | | | | | |
| | | | | | | |
| Noto: | All SKC tubes ar | e manufactured | with a space hefo | re the enring a | lace wool primary co | action congrator |

Note: All SKC tubes are manufactured with a space before the spring, glass wool, primary section, separator, backup section and final separator.



| В. | Separators (Retaining plugs): Indicate which type for each separator required. | | | | | | | |
|-------------------|--|---|--------------------------|--|--|--|--|--|
| | Separator | Type Required | Separator | Type Required | | | | |
| | Α | Select One | С | Select One | | | | |
| | В | Select One | D | Select One | | | | |
| Stainle | ss Steel Spring F | Retainer Required? Select One | | | | | | |
| C. Sorl | bent: | | | | | | | |
| Note: I | f sorbent will be | e supplied by customer, SDS may | be required. | | | | | |
| Section Sorben | n 1: t type Select One | | Bed weight | | | | | |
| Section Sorben | n 2: It type Select One | | Bed weight | | | | | |
| Section Sorben | n 3: t type Select One | | Bed weight | | | | | |
| Sorben | t preparation or t | treatment: | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| | | | | | | | | |
| 6. Purd | chasing Terms a | and Authorization: | | | | | | |
| Ехре | edited service may | y be available for a fee. Contact Custom | er Support to make arra | e for any custom media is 4-6 weeks ARO. ngements. Sorbent bed weights are \pm 10%. ts call 800-752-8472 (1-724-941-9701). | | | | |
| not b | Terms: Purchase orders for solvent desorption sorbent tubes must be accompanied by this form. Purchase orders manned be changed or cancelled after they are received by SKC Inc. Custom order solvent desorption sorbent tubes are not returnable. See the SKC Limited Warranty and Return Policy at http://www.skcinc.com/warranty. | | | | | | | |
| | SKC reserves the right to cancel this order if it is determined that SKC cannot produce a quality solvent desorption sorber tube according to the specifications provided. | | | | | | | |
| | | signature of the undersigned confirmes to the custom order terms. | s that the specification | s supplied on form 9560, Revision 2212 | | | | |
| Name (| (please print or ty | /pe) | | | | | | |
| Signatu | ure | | Date | | | | | |
| Compa | ny Name | | | | | | | |
| Contac | t Telephone Nun | nher | Fmail | | | | | |